



PROFESSIONAL QUALITY, FOCUS ON VARIOUS
OXYGEN CONCENTRATORS AND OXYGEN RELATED PRODUCTS



oxygen care
specialists

Imported & Distributed by:



SOUTHSIDE BIOMEDICAL SALES AND SERVICES CO. LTD.

 #903 9TH FLOOR PHILFLEX BAY CENTER,
CORAL WAY DRIVE, PASAY CITY

 # 8-461-3468/ 0917 728 0229

 info@southsidebiomedical.com





COMPANY PROFILE




Longfian Scitech Co.,Ltd, located in 2F&3F, East Section, Building 12, Power valley pioneer park, NO.369 Huiyang street , 071051,Baoding, Hebei Province, China, is a professional manufacturer, specialized in the development, production and export of medical serial oxygen concentrators.In 2016,it established second factory named"Dongguan Medjoy Medical Technology Co.,Ltd",located at Changping Town,Dongguan,Guangdong ,China.The company is possessed with its own solid R&D team, who has the independent developing ability of the hard&software unit of the product and won the matured technology of PSA application. Presently the company has succeeded in the production of serial oxygen concentrators from 1LPM~60LPM, which have been widely applied in homecare, hospitals, clinics, beauty salon, sports, industry and high altitude area etc.

In the meantime, some oxygen related products, such as ozone, oxygen cocktail device, are developed also upon the demands of the market.

In order to meet the international quality requirement, some products have got FDA certificate and TUV approval with CE certificate and exported to United States, France, Spain, Germany, Russia, India, Singapore, Turkey, Thailand, Colombia, Mexico, Argentina, Kenya, Ethiopia, Nigeria etc. more than 160 countries and regions, enjoy a high reputation in international market.

Our goal is professional in quality, focus on various oxygen concentrator & oxygen related products, and continue to bring you not only reliable, quality products but also bring you a better life.

CERTIFICATE



HEMOCARE OXYGEN CONCENTRATOR

MODEL: JAY-2/JAY-3



Pediatric
(neonatal)
machine



Small & light



0-2 L flow meter
Minimum Scale 0.125L/Min



Application

- JAY-2: The flow rate is 0-2L. The minimum scale can be 0.125L/min. It is designed specially for infants.
- JAY-3: The flow rate is 0-3L. For home care and health care

Specification

Flow Rate	0-2LPM / 0-3LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	≤45db
Power Consumption	≤200VA / ≤250VA
LCD Display	Present working time, accumulated working time
Alarm	Power Failure alarm, Low purity alarm (Optional)
Net Weight	11Kgs / 12Kgs
Size	435×240×410mm

HEMOCARE OXYGEN CONCENTRATOR

MODEL: JAY-3AW/JAY-5AW



Low noise,
takes up less space



USA imported humidifier,
6PSI pressure release



Innovative big
LCD display

Optional

- Nebulizer (Atomization)
- Low Purity alarm
- Low purity alarm with low flow alarm
- High temperature alarm
- Pulse oximeter
- Maintenance alarm
- Reset Switch
- Remote control monitoring system (WIFI/Bluetooth)

Specification

Flow Rate	0-3LPM / 0-5LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	≤40db
Power Consumption	≤300w / ≤350w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	16Kgs / 16 Kgs
Size	348×280×510 mm

HEMOCARE OXYGEN CONCENTRATOR

MODEL: JAY-3BW/JAY-5BW



Higher outlet pressure to ensure proper flow delivery



Low maintenance, easy to care for



Optional pulse oximeter

Optional

1. Nebulizer(Atomization)
2. Low Purity alarm
3. Low purity alarm with low flow alarm
4. High temperature alarm
5. Pulse oximeter
6. Maintenance alarm
7. Reset Switch
8. Remote control monitoring system (WIFI/Bluetooth)

Specification

Flow Rate	0-3LPM / 0-5LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	≤43db
Power Consumption	≤300w / ≤350w (≤540w)
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	16Kgs / 18Kgs
Size	360×300×600mm

MEDICAL OXYGEN CONCENTRATOR


MODEL: JAY-5/JAY-8/JAY-10



DUAL FLOWS
USA imported humidifier
6PSI pressure release



Innovative "cooling" technology to improve reliability and life of the unit



Big LCD display
24H*365D work continuously

Optional

1. Nebulizer(Atomization)
2. Low Purity alarm
3. Low purity alarm with low flow alarm
4. High temperature alarm
5. Pulse oximeter
6. Maintenance alarm
7. Reset Switch
8. Remote control monitoring system (WIFI/Bluetooth)

Specification

Flow Rate	0-5LPM / 0-8LPM / 0-10LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	≤50db
Power Consumption	≤550w / ≤550w / ≤880w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	26 Kgs / 26 Kgs / 27 Kgs
Size	365×375×600mm